

WAS[®]

MAKING VEHICLES SPECIAL



Comfortable and Ergonomic:
**The New MAN TGE
WAS 300 Ambulance.**

WAS | 300 AMBULANCE 4.25 T

The New MAN TGE WAS 300, 4.25 T: A Efficient and Well-Organised Vehicle for Any Job.



SHORT DESCRIPTION

This ambulance impresses with its smoothly finished surfaces and tidy interior, which allows for optimal use of space. Medical equipment is safely stored in cabinets to ensure the highest level of hygiene. The patient is the focus, while the safety of the paramedic is always guaranteed. Well-thought-out storage space in the right places, which does not interfere with the paramedic's movements, ensures comfort and enables efficient work. The result is a pleasant workplace that offers everything you need. The MAN TGE base vehicle impresses with its robust construction and high level of comfort. The driver and paramedic benefit from a spacious workplace with high-quality seats and personalised leather covers. A strong, economical engine and guaranteed durability make this vehicle the ideal choice.



Ergonomically designed storage space for the safe storage of a carrying chair, scoop stretcher and oxygen cylinders.



The sliding door on the right-hand side offers spacious access and features a closing aid. The attached foil on the integrated window ensures that the patient's privacy is always protected. Full-length sidebars provide safe access while effectively protecting the body from damage.

EXTERIOR

All-round LED work lights
 Integrated multi-colour warning system at the front and rear, type DON-light
 LED front flashers with silicone lens (blue)
 LED mirror flashers with silicone lens (blue)
 LED mudguard flashers with silicone lens (blue)
 LED flashers with silicone lens on rear doors (blue)
 Acoustic warning system integrated in the radiator grille
 LED entrance lighting
 Rear view camera and surveillance sensors integrated in the rear spoiler
 Sidebars with integrated step
 OOV decals according to SOVAM regulations

INTERIOR

Stryker Powerload stretcher system
 Various storage compartments on the partition wall:
 full-height cupboard with pull-out drawer, cupboard with tablet and ECG drawer, integrated ampoule drawer
 Attendant seat on the right, adjustable longitudinally
 Attendant seat on the partition wall, seat can be folded up
 Multifunctional ceiling-high cabinet at the rear right for vacuum mattress and Lucas 3
 Multifunctional ceiling-high cabinet at the rear right for vacuum mattress and Lucas 3
 Oxygen supply with various oxygen sockets
 230 V inverter
 230 V and 12 V sockets
 CAN bus central electrics 12 V
 Heat exchanger and climate control, controlled via central touch display
 Thermobox with temperature display
 Roof storage compartment with LED lighting in partition wall
 External storage compartment behind the left sliding door for stowing: scoop stretcher, carrying chair, 2 x 10 l oxygen
 Left-hand hanging cupboard with 3 flaps and variable compartment subdivision
 Left-hand storage cupboard with 3 flaps, to be opened upwards
 LED interior lighting with emergency lighting and blue trauma lighting
 LED doctor's spotlight, switched via control unit
 Electric heater, controlled via CAN bus system
 Roof vent for ventilation, 12 V
 Interior camera

DRIVER'S COMPARTMENT

Navigation and patient data shown on original MAN display
 Second display for navigation and patient data on the co-driver's side
 Radio equipment (C2000)
 Foot switch for warning system
 Centre console between driver and co-driver's seat with integrated cool box
 Touch display positioned for easy access
 Central electrics 12 V, installed behind the driver's seat
 Charging technology and inverter, installed behind the co-driver's seat
 Comfortable seats: with seat cushion extension, leather cover including logo and armrests



Multifunctional cabinet with spacious and secure storage that allows for quick and easy access.



High-quality and modern cockpit with an unobstructed view of the road.

The MAN TGE WAS 300 Ambulance 4.25 T: State of the Art in Ergonomics and Comfort.



Tidy and easy-to-disinfect patient compartment with a wide access on the right side.

SPECIAL FEATURES

Multifunctional storage compartments, accessible from the inside and outside

Comfortable and versatile companion seats

Automatic switchover unit for oxygen supply

Heating and air conditioning in the patient compartment

BASE VEHICLE

Model	MAN TGE
Vehicle Type	3.180 L3H3, 137 kw (180 hp), 4 x 2
Transmission	8-speed automatic transmission
Equipment	Storage package: storage system in the roof with two 1 DIN slots and reading light LED headlights 'Climatronic' air conditioning MAN Media VAN package Business navigation Rear doors with increased door opening angle Driver and passenger seats 'Comfort Plus' Engine idle-run switch Closing aid, sliding door on the right and left Second battery with high-current cut-off relay and battery monitoring
Wheel base	3640 mm
Maximum total weight	4250 kg



Optimal seating position close to the patient, in a safe environment.



Simple and robust design and easily accessible storage space.

The concept is transferable to other vehicle types. The equipment is an example and can be customised.



Wietmarscher Ambulanz- und Sonderfahrzeug GmbH | Darwinstr. 11 | 48488 Emsbüren, Germany
Phone +49 5903 93201-100 | Fax +49 5903 93201-602 | info@was-vehicles.com | www.was-vehicles.com