

**WAS**<sup>®</sup>

MAKING VEHICLES SPECIAL



The New Swedish Standard:  
**VW Crafter WAS 500 Euro-Lans.**

WAS | 500 EURO-LANS

# The New VW Crafter WAS 500 Euro-Lans: Modern Interior Design with Very Compact Dimensions.



## SHORT DESCRIPTION

The VW Crafter WAS 500 Euro-Lans ambulance presents all modern standards of the WAS brand combined with a highly ergonomic interior. The patient compartment offers best sanitary working conditions and enough space for a comfortable treatment of the patient. All required materials are not only stored functionally but are also within reach. Safety of emergency services personnel has the highest priority therefore the defibrillator and the main emergency case are mounted on a sled which can be pulled in direction of the paramedic. In this way the paramedic has all his attention for the patient without being forced to leave his seat. The ergonomic design with enclosed streamlined panels between cab and box body ensures optimum aerodynamics, thus optimising the vehicle's fuel consumption.





## EXTERIOR

LED working lights  
 Special front and rear LED warning system, with Standby W3 bar  
 LED front flashers in white and blue, type Standby  
 LED flashers integrated in the side view mirrors  
 LED flashers integrated in the rear bumper  
 Warning system FSC 100 Composite, mounted behind grill  
 LED flasher blue/amber at the top of box body  
 LED flasher amber integrated in rear warning system the left and right side wall  
 Integrated rear LED warning system  
 Rear LED working lights  
 360° camera system, type Brigade  
 Roof spoiler above driver's cab  
 Fold-out side step, automatically triggered by entry door (mechanical powered)  
 Outside step at rear doors  
 Outside walls made of polyester-coated aluminium  
 LED entry lighting in all entrance areas  
 LED vehicle taillamps, side marker lights and clearance lights



## INTERIOR

Attendant seat on the partition wall, with seat heating and seat recognition  
 Storage for communication devices, equipped with an extendable tray  
 Special slidable carrying system for defibrillator and emergency bag  
 LED interior lighting with emergency light and blue trauma light  
 Two attendant seats with integrated seat heating and seat recognition, mounted on fully coated wheel houses  
 Roof storage compartment above bulkhead window, with integrated LED lighting  
 2 air outlets at the left front cabinet  
 4 oxygen outlets at the left front cabinet  
 Upper cabinet with transparent flaps  
 Working surface with integrated keyboard pull-out table  
 State of the art (Standby) control panel  
 Compartment for oxygen/kallinox cylinders, with additional storage above  
 Document holder positioned within reach of paramedic (left side)  
 Rail mounting system integrated in left side wall  
 Ceiling center with integrated infusion box, oxygen/kallinox outlets and radio speaker  
 Compatible with all current European stretcher systems  
 Automatic climate control for AC system and auxiliary heating  
 Oxygen outlets  
 12 V outlets  
 230 V outlets, with power supply during transport



# The New VW Crafter WAS 500 Euro-Lans: Setting New Standards in Safety and Ergonomics.

## DRIVER'S COMPARTMENT

Centre console between driver's and co-driver's seat housing equipment and A4 documents  
 Large control panel with touch screen  
 Monitor for 360° camera system  
 Loading technology and 230 V inverter easy to access behind co-driver's seat  
 12 V CanBus system and 230 V central electric behind driver's seat  
 Two Maglite charging cradles at front of centre console  
 Two universal helmet holders in the roof storage

## BOX BODY

Model                   Emergency ambulance type B or C  
 Length: 3490 mm  
 Interior  
 Width: 1960 mm  
 Height: 1800 mm

Doors

- Side hinged door on the right
- Rear double hinged doors
- All doors and compartments with central locking

Compartments

- Large access compartment door, front right side
- Rear access compartment, right side, for large equipment like scoop stretcher etc.
- Compartment front left side for LUKAS or breathing bag

Vacuum mattress, fixation, carry chair, scoop stretcher, mini stretcher and bags are fixed and can easily be removed from vehicle

Compact and safe: oxygen-, air- and gas bottles inside

## SPECIAL FEATURES

Hi-Tech CanBus technology  
 360° camera system  
 Attendant seats with heating and seat recognition  
 Slidable carrying system for defibrillator and emergency bag  
 230 V inverter with 1800 Watt  
 Fold-out side step, automatically triggered by entry door (mechanical powered)  
 Airsuspension for ergonomic loading of the patient

The concept is transferable to other base vehicle types.

The equipment is an example and can be customised.

## BASE VEHICLE

Model	Volkswagen Crafter 35
Vehicle type	2.0 TDI EU 6 BlueMotion
Engine power	130 kW/177 hp
Gear box	Automatic transmission
Wheel base	3640 mm
Maximum GVW	4300 kg

